
BONE MANAGEMENT® CATALOGUE

BONE AUGMENTATION TECHNIQUES

– YOUR KEY
TO SUCCESS

EDITION III

GLOBAL BONE MANAGEMENT® & RESTORATIVE NETWORK

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When developing our surgical systems and instruments, we are not only concerned with the highest quality of "Made in Germany", but also with the correct application of these systems and instruments. Since the professional preparation of the bony implant site is essential for successful implant placement and management, MEISINGER has developed a multi-level training concept in cooperation with leading scientists and practitioners. With the Bone Management® Academy, both experienced users and beginners have the opportunity to deal intensively with the topics of bone augmentation, implantology, content related to practice management as well as providing the chance to further their education in this field.

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Our main motivation is on the continuously increasing added value for participants. The focus is therefore not on the systems developed by MEISINGER, but on professional, practical, and patient-oriented techniques that can be learned in the courses. High-quality and precisely manufactured products play an essential role in positive and successful patient outcomes. However, the correct and proficient application of these tools is the famous "key to success" and thus also the best way to a satisfied patient.

That is why we do not shy away from our competition, but deliberately involve other dental companies in order to enable our participants to form a balanced perspective, and to discuss and compare competing treatment methods. If the participants find certain techniques to be good, the products can firstly be intensively tested and practiced in the workshops. Additional work shadowing and supervision offer further opportunities to consolidate what has been learned and thus make better informed decisions.

The Young dental experts camps, which take place several times a year, also offer orientation aid for students, interns, and residents in their future specialization choice. Aspiring surgeons can also acquire solid basic knowledge in the various courses and workshops. The right suturing and cutting techniques are just as important as the right understanding of soft and hard tissue. Various augmentation techniques can also be learned and practiced within the curriculum. The Bone Management® Academy also includes complex techniques such as the Carrot and Shell techniques developed by Prof. Dr. Fouad Khoury or the Sausage technique according to Dr. Istvan Urban.

MEISINGER USA Headquarters, Colorado

MEISINGER Germany Headquarters, Neuss

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22 BONE MANAGEMENT® MASTER-LINE

22	BMCBA	Master-Core Basic
23	BMCPR	Master-Core Professional
24-25	BKM00	Master-Mill
26	BMPBA	Master-Pin-Control Basic
27	BMP00	Master-Pin-Control
28-29	BMPPL	Master-Pin-Control Plus
30	BMPXB	Master-Pin-X-Control Basic
31	BMPXU	Master-Pin-X-Control Upgrade
32-33	HF020	Master Surgery Kit

38 BONE MANAGEMENT® KHOURY-LINE

38	BMSTP	Soft Tissue Punch System
40	BIR00	Implant Removal Kit
41	BIE00	Implant Extraction Kit
42-43	BTE00	Trephine Ejection Kit
44	BMSTB	Micro Screw T5 System Basic
45	BMST0	Micro Screw T5 System
46-47	BMTSP	Micro Screw T5 System Plus
48	BMSBA	Micro Screw System Basic
49	BMS00	Micro Screw System
50-51	BOST0	Ost-Tray
52	BSIN0	Sinus-Tray
53	BSOT0	Soft-Tissue-Tray
54	915DM	Diamond Disc
54	943DC	Diamond Disc
54	M1159	Bone Block Tweezers
55	BS005	Smartscraper
55	MBHD1	Micross
56	BS003	Safescraper® TWIST straight Version
56	BS004	Safescraper® TWIST Curve angled Version

58 GBR, BONE TRANSFER & OSTEOSYNTHESIS

58	BTS00	Tenting Screw System
59	USL600	Umbrella Screw
60	GBR00	My First GBR Kit
62-63	CTR00	Transfer-Control
64	CTR10	Transfer-Ring-Control I
65	BBR00	Transfer-Ring-Control II
66-67	CTRPL	Transfer-Control Plus
68	BTX00	Screw System TX
69	BTXPR	Screw System TX Professional
70	7120	Trephine Basic Kit
71	7121	Trephine Kit
72	7122	Trephine Kit Long
73	7140	Punch Basic Kit
73	7150	Saw Basic Kit

74 SINUSLIFT

74	CCLBA	Crestal-Lift-Control Basic
75	BCL00	Crestal-Lift-Control
76	BTL00	External-Lift-Control
77	CSL00	Internal-Lift-Control
78-79	BLIPL	Lift-Control Plus
80	BSKSL	Surgical Kit 1
81	BSK02	Surgical Kit 2

82 IMPLANT SITE PREPARATION

82-83	BIB00	Implant Basic Kit
84	BSK03	Surgical Kit 3
84	BSK04	Surgical Kit 4
85	ALV18	Alveoplasty Surgical Kit
86-87	BGS01	Guided-Drill-Stop-Control

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Note: Not all products contained in this catalogue are available in all countries.

Benex®-Control

Root Extraction System |
developed with Dr. Benno Syfrig

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Benex®-Control is outstandingly suitable for the controlled and safe extraction of dental roots and dental root fragments, so the conventional extraction method using forceps and leverage can be avoided. Thanks to the innovative construction of the extractor, the root can be easily removed in a particularly controlled manner without damaging bone or soft tissue. To guarantee an optimal vertical transmission of force, there must be a flat area providing a supporting surface to the extractor. The quadrant support was developed specifically to provide this optimal support while also protecting existing dentition and providing support for the extractor in case of edentulous circumstances.

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
Fig. BE001 015
gerade | straight

¹ Teflon® is a registered trademark of DuPont, USA

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AT A GLANCE

- Safe and controlled extraction of roots
- Possibility of immediate implant placement
- Maximum preservation of the vestibular bone wall due to vertical extraction
- Optimal protection of bone and soft tissues
- Applicable to small root fragments
- Osteotomy can be prevented
- High patient acceptance

Developed with
Dr. Benno Syfrig
 Luzern, Switzerland


Get 'em all out Kit NEW

Root Extraction System | according to
Dr. Felix Roth & Dr. Fabian Schiml

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The Get 'em all out Kit consists of two surgical burs (two of each) and five hand instruments and has been specially designed for minimally invasive and predictable tooth removal in simple and complex cases. Thanks to the coordinated procedure, it is suitable for both beginners and advanced users. Using the hand instruments and burs included, crowns can be removed, roots can be divided and even root tips can be exposed by intra-alveolar osteotomy. The focus is on maximum bone preservation with regard to subsequent implant placement. The ergonomic handle shape of the luxation instruments ensures controlled and safe handling. Residual fragments can be removed easily and safely using the fragment ejector and tweezers. A curette is used for final cleaning of the extraction socket.

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AT A GLANCE

- Sophisticated assortment of instruments for effective and structured extraction (for both beginners and advanced users)
- Safe handling thanks to the ergonomic handle shape of the luxation instruments
- High flexibility thanks to two different luxation instrument variants (straight, angled)

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RTRO Kit **NEW**

Root Tip Resection and Osteotomy |
developed with PD Dr. Dr. Florian Bauer

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The RTRO Kit by PD Dr. Dr. Florian Bauer contains a sophisticated selection of diamond instruments and carbide burs for performing root tip resection (apicoectomy and osteotomy). The HM141A and 801 instruments are used to first expose the root apex by removing the overlying bone tissue. Whereas the 837L, 850G, HM165 and HM166A instruments are applied in the subsequent removal of the root tip as well as of inflamed tissue in its vicinity.

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AT A GLANCE

- Sophisticated combination of high-quality diamond and carbide burs for apicoectomy and osteotomy
- All rotary instruments for exposure and tissue removal in one single kit
- Safe and efficient exposure and removal of the root tip by high-precision instruments

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The Root Elevator Kit Samarani consists of 6 elevators and a small hammer. It was specially developed by Dr. Rawad Samarani for the safe and novel POP-OUT root extraction and is used in combination with Kit BRSS2 for root sectioning. The root fragments are removed from the socket by light, horizontal hammering using the newly designed elevators. The load is applied to the fragment by means of the sophisticated angulation of the elevator's working part, which is specifically adapted to the respective anatomical position where the extraction takes place. The laser markings on the elevators indicate the intraoral position of application (buccal, palatal / lingual). Additional elevators are included for special cases where the root to be extracted is below the level of the alveolar bone. The instruments can be conveniently sterilized directly in the tray supplied and transported safely. The elevators do not put loads on the adjacent tissue as it is the case with conventional extraction techniques. The innovative load application significantly reduces the risk of injury of soft tissue, bone or adjacent teeth. The high level of safety of this extraction technique is a decisive advantage, particularly when immediate implant placement is intended.

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The Root Sectioning Kit Samarani comprises 6 selected tungsten carbide instruments for sectioning multi-rooted teeth into fragments. By combining this kit with BRSS1, extraction of roots can be performed as follows:

- If the tooth is multi-rooted, sectioning of the tooth into separate root fragment by using burs from BRSS2
- Selection of an elevator from BRSS1 according to the anatomical position of the root
- Intraoral positioning of the elevator at the root
- Atraumatic mobilization of the root by gently hammering on the end of the elevator

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AT A GLANCE

- Straightforward and atraumatic root extraction and thus preservation of soft tissue and alveolar bone
- Specifically shaped and angled elevators for innovative load transmission leaving the surrounding tissues untouched
- Small hammer for maximized patient comfort
- Preservation of extraction socket as a prerequisite for immediate implant placement and satisfactory esthetic results

Hartmetallinstrumente Tungsten carbide burs						
Fig.	HM166RX*	HM166RX*	HM141F*	HM151*	HM151*	HM1*
Shank ¹	205	206	206	316	317	316
Size ²	021	021	010	016	016	008
Length mm	11.0	11.0	0.88	10.8	10.8	0.71
* mm	2.1	2.1	1.0	1.6	1.6	0.8

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



Crest-Control is a Horizontal Bone Splitting System which allows for a predetermined and minimally invasive widening of the horizontally resorbed alveolar ridge particularly in the distal lower jaw. With the well-proven Horizontal Spreaders, the alveolar ridge can be expanded up to a maximum of 5 mm, allowing an implant of every common system to be inserted in the widened segment.

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


AT A GLANCE

- Controlled widening of horizontally resorbed bone
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all common implant systems
- Minimally invasive treatment concept
- Gentle widening of the alveolar ridge

Split-Control 12 mm

Bone Condensing and Bone Spreading System

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





Split-Control 12 mm is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

Split-Control 12 mm provides an extensive range of instruments for a controlled dilatation with spreaders with a length of 12 mm.

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



AT A GLANCE

- Controlled spreading of horizontally resorbed bone
- Gentle bone condensing with the aid of non-cutting Archimedes screws
- Increased primary stability of the inserted implants
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all common implant systems

Split-Control Plus

Combination of Split-Control and Crest-Control

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Split-Control Plus is a combination of the proved Split-Control and Crest-Control systems, so it offers an extensive range of instruments for controlled bone spreading and condensing. It contains both specially designed screw-like condensing and spreading instruments (Spreaders) and horizontal spreaders that allow for controlled and standardized spreading of horizontally resorbed bone. Due to the special geometry of the screw-like spreaders, cancellous bone is gently condensed while the Horizontal Spreaders allow for a dilatation up to a width of 5 mm, so the bone is optimally prepared for the insertion of every common implant.






Anwendungsvideo
Application Video

Art.-No. CSPPL



© Dr. Boyer & Drs. Kötter, © Dr. Raymond Kimsey, D.D.S.

	Vorkörner Initial bur	Diamantierte Sagen Diamond coated saws		Osteotomieinstrumente Osteotomy instruments		Pilotbohrer Pilot burs		Erweiterer 12 mm Expansion burs 12 mm	
									
Fig.	186RF	231DC*	231DC*	HM331L**	HM254E**	A1001	D1001	A2003	D2003
Shank ¹	204	204	204	205	205	206	206	205	205
Size ²	018	100	130	010	012	010	018	023	030
Length mm	12.0	0.3	0.3	5.5	6.0	18.0	18.0	12.0	12.0
 mm	-	-	-	-	-	-	-	1.30	2.03
 mm	1.8	10.0	13.0	1.0	1.2	1.0	1.8	2.3	3.0

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AT A GLANCE

- Combined System of Crest- and Split-Control
- Possibility of combining bone splitting and spreading techniques
- Horizontal spreaders for gentle widening of the alveolar ridge by up to 5 mm
- Non-cutting Archimedes screws for gentle bone condensing
- Preparation of the bone for insertion of all common implant systems

BONE MANAGEMENT® MASTER-LINE

Developed with
Dr. Istvan Urban, DMD, MD, PHD

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MEMBER OF THE GLOBAL
BONE MANAGEMENT® & RESTORATIVE NETWORK

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The Master-Core Basic is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trepanes with different lengths and diameters. Due to the different lengths, the trepanes have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trepanes' working parts enable additional depth control during the extraction. The black coating of the trepanes allows for glare-free work and the optimal visibility of the depth markings.

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AT A GLANCE

- Simple and safe extraction of bone cylinders
- Optimal protection of anatomical structures thanks to depth-limited trepanes
- Depth markings on the trepanes for additional depth control
- Extraction of bone cylinders of different diameters and lengths
- Black coating on the trepanes for glare-free work

	Vorkörner Trepane Initial bur Trepanes		Bohrerverlängerung Drill Extension
Fig.	DV229	DV229	BV025
Shank ¹	204	204	205
Size ²	050	070	-
Length mm	2.0	2.0	25.5
mm	5.0	7.0	-
mm	5.7	7.7	-

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
The Master-Core Professional is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trepanes with different lengths and diameters. Due to the different lengths, the trepanes have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trepanes' working parts enable additional depth control during the extraction. The black coating of the trepanes allows for glare-free work and the optimal visibility of the depth markings. To be able to respond more precisely in specific indications, the Master-Core Professional contains a broad range of instruments.

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Master-Mill

Surgical Bone Mill – used by Prof. Istvan Urban

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The Master-Mill is a surgical bone mill used for crushing autologous bone blocks and thus preserving small autologous bone chips. With the help of these bone chips, bone defects of all shapes can be optimally filled. The individual components of the Master-Mill are made of high quality, surgical stainless steel. All parts feature an optimal fit and a stable geometry, allowing the Master-Mill to be easy to handle with ease of assembly and intuitive use.

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AT A GLANCE

- Crushing of autologous bone blocks
- Augmentation of bone defects of all shapes
- Individual components made of high quality, surgical stainless steel
- Stable geometry
- Easy and intuitive handling

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
Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5-ELI (Titanium Grade 23) and have a groove which increases the surface of the pins. Due to the groove and a unique head-design, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

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AT A GLANCE

- No bending, no breakage - Safe Fixation of Membranes with Pins made of Titanium Alloy Grade 5-ELI (Titan Grade 23)
- Sharp tip and high stability of the pins allow for precise placement of the pins also in cortical bone and at an angle
- Easy Removal from the tray and placement of the pins due to the specially developed Pin-Holder
- Instruments for the placement of Bleeding Points for optimal healing conditions

	Vorkörner Initial burs (2 Stück. / 2 pcs.)	Spiralbohrer Twist drills (je 2 Stück. / each 2 pcs.)	Pins (10 Stück. / 10 pcs.)	
Fig.	186RF	203RF	203RF	MP10*
Shank ¹	204	205	205	-
Size ²	018	006	008	-
Length mm	12.0	7.0	7.0	3.65
mm	-	-	-	0.87
mm	1.8	0.6	0.8	0.95
Head-Ø mm	-	-	-	2.51

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Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The Master-Pin-Control Plus contains two different types of pins. While one pin has a particularly flat head design and is therefore, for example, ideal for use under thin soft tissue, the other pin has a larger head design that makes it particularly easy to find the pins under the soft tissue after successful healing. All included pins consist of titanium alloy Grade 5-ELI (Titan Grade 23) and have a groove which increases the surface of the pins. Due to the groove, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone. For the simple and safe placement of bleeding points in the cortical bone, the Master-Pin-Control Plus system also contains the twist drill 203S, which has an automatic depth stop.

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
AT A GLANCE

- No bending, no breakage - Safe Fixation of Membranes with Pins made of Titanium Alloy Grade 5-ELI (Titan Grade 23)
- Two different types of pins for optimal care in every soft tissue situation
- Sharp tip and high stability of the pins allow for precise placement of the pins also in cortical bone and at an angle
- Easy Removal from the tray and placement of the pins due to the specially developed Pin-Holder
- Twist drill 203S with automatic depth stop for the simple and safe placement of Bleeding Points for optimal healing conditions

Master-Pin-X-Control Basic NEW

Pin System for Membrane Fixation | developed with Dr. Istvan Urban

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The Master-Pin-X-Control Basic (BMPXB) comprises ten pins with thread (MP30), which are contained in a pin box. The matching screwdriver (TL0T1) for manual insertion is also included. The system also includes two twist drills with depth stop (203S) and two different drill bits for pre-drilling (203RF).

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AT A GLANCE

- Secure fixation of membranes thanks to new pins with thread made of biocompatible titanium alloy Grade 5-ELI (Titan Grade 23)
- Easy and safe transfer and optimum force transmission during insertion thanks to Torx interface
- Drill with depth stop for safe placement of bleeding points for optimised bone regeneration
- Secure attachment of membranes now also in hard-to-reach areas (lingual)
- Thank to self-tapping thread generally no need to pre-drill (only recommended for high bone hardness)

	Spiralbohrer Twist drills (je 1 Stck. / each 1 pcs.)		Spiralbohrer mit Stop zum Dekortikalisieren Twist Drill with Stop for decoricating (je 2 Stck. / each 2 pcs.)	Pins Pins (10 Stck. / 10 pcs.)
Fig.	203RF	203RF	203S	MP30*
Shank ¹	205	205	204	-
Size ²	006	008	012	-
Length mm	7.0	7.0	4.0	3.2
mm	0.6	0.8	1.2	-

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The Master-Pin-X-Control Upgrade (BMPXU) is a Master-Pin-X-Control system reduced to the minimum required components, which is the most compact option for the user. The pin box contains 20x MP30 pins. A TL0T1 screwdriver is also supplied. Note: The pin box of the MP10 pin without thread is not compatible with the MP30 pins!

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AT A GLANCE

- Secure fixation of membranes thanks to new pins with thread made of biocompatible titanium alloy Grade 5-ELI (Titan Grade 23)
- Easy and safe transfer and optimum force transmission during insertion thanks to Torx interface
- Secure attachment of membranes now also in hard-to-reach areas (lingual)
- Thank to self-tapping thread generally no need to pre-drill (only recommended for high bone hardness)

Pins Pins (20 Stück / 20 pcs.)	
	
Fig.	MP30*
Length mm	3.2

Schraubendreher manuell Screwdriver manual illustrated 1:2	
	
Fig.	TL0T1
Length mm	148.0

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This surgical kit contains 18 hand instruments, 2 mouth props and a drill stand with six slots for rotary instruments. Developed by Dr. Istvan Urban, it is optimally suited to be used especially in minimally invasive interventions for hard and soft tissue regeneration. Within the tray supplied, the instruments can be sterilized and transported safely. Instruments included with a black finish show a significantly reduced reflection of light.

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AT A GLANCE

- Sophisticated product combination in one kit for performing minimally invasive surgical interventions according to Dr. Istvan Urban
- Sterilization and transport within the supplied tray
- Black finish on selected instruments for reduced reflection of light

THE ART AND SCIENCE OF
TISSUE THERAPIES

MASTER CLASS ON TOUR - 2026

DR.
ISTVAN URBAN

PROF. DR.
FOUAD KHOURY

29. / 30.05.2026
PARIS | FRANCE

12. / 13.06.2026
MILANO | ITALY


SAVE THE
DATES

04. / 05.09.2026
DÜSSELDORF | GERMANY

06. / 07.11.2026
WARSAW | POLAND

MEET US AT EAO **MONACO**

THE IMPACT OF TIME IN IMPLANT DENTISTRY

Speaker:
Prof. Dr. Fouad Khoury, Germany

Title:
**CURRENT CONCEPTS FOR SIMPLE
AND SECURE BONE AUGMENTATION**

Abstract:

This hands-on course will cover different aspect of bone grafting procedures with intra oral harvested bone as well as the management of the soft tissue. Special focus will be on a minimal invasive but very effective technique, the carrot technique, that is predictable, reproducible and easy to implement in clinical practice. This approach reduces morbidity, shortens treatment time, and eliminates the need for additional donor sites or foreign materials. Step-by-step guidance will be provided to enhance long-term success in implantology through biologically optimized grafting procedures.

Vita:

- Born in Lebanon
- 1978 DMD, St. Joseph University, Beirut
- 1978-1979 Department of Oral and Maxillofacial Surgery of the University of Freiburg, Germany
- 1979-1994 Assistant Professor and later Associate Professor at the Dep. of Oral and Maxillofacial Surgery of the University of Muenster, Germany
- 1984 Diploma in Oral Surgery and Doctorate in Medical Dentistry
- 1988 Habilitation and senior surgeon at the dep. of Oral and Maxillofacial Surgery of the University of Muenster
- 1994 Professor at the University of Muenster

**MEISINGER
2X WORKSHOPS
AT EAO:**

BONE MANAGEMENT® KHOURY-LINE

Developed with
Prof. Dr. Fouad Khoury

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Micro Screw T5 System / Micro Screw System

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
Soft Tissue Punch System

NEW

For Biopsies or tissue grafts |
developed with Prof. Dr. Fouad Khoury

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The Soft Tissue Punch System contains mucosal punches in five different diameters (5-13 mm). The 9 mm, 11 mm and 13 mm sizes have an extended shaft in order to facilitate the access to areas deeper in the oral cavity. The punches are used to obtain biopsies or soft tissue grafts, for example for immediate implant placement or in order to create access to a previously placed and healed implant. Thanks to precise manufacturing and high sharpness, the mucosal punches have excellent cutting properties and thus reduce the risk of tissue trauma.

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
AT A GLANCE

- Selection of high-quality mucosal punches for the harvesting of soft tissue grafts
- Extended shank length to improve access to more remote areas in the oral cavity
- High sharpness and thus a reduced risk of tissue trauma

myplant

T W O

TWO PERFECT SOLUTIONS!

OKTAGON®

MEISINGER
IMPLANTS

Meisinger

Implant Removal Kit NEW

developed with Prof. Dr. Fouad Khoury

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The Implant Removal Kit is designed for the simple and safe removal of implants from the jaw. The threaded pins included in the kit contain two threads and can be screwed with one side into the internal thread of the implant to be removed. With the help of the removal shanks, which are screwed onto the other side of the threaded pins, the implant and the removal shank are locked together, allowing the torque to be transferred from the ratchet to the implant. After loosening, the implant can be unscrewed without damaging the tissue. After thorough medical assessment of the implant bed, a new implant with an identical diameter might be inserted.

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The Implant Extraction Kit allows the guided extraction of osseointegrated implants. Thanks to the innovative yet simple guiding technique, the trephine diameter can be adjusted as closely as possible to the implant diameter. Thus, the extraction can be performed in a way that preserves most of the surrounding bone tissue. The kit contains five different sleeves, five different pins with threads and ten trephines, which differ in both diameter and length. For extraction of an osseointegrated implant, the matching threaded pin is screwed into the internal thread of the implant. The trephine is then selected based on the length and diameter of the implant. The matching sleeve is then attached to the threaded pin on top of the implant and serves as a guide for the trephine during extraction.

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Trephine Ejection Kit


Trephine System for Extraction of Cylindrical Bone Grafts | developed with Prof. Dr. Fouad Khoury

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
The Trephine Ejection Kit was developed together with Prof. Dr. Fouad Khoury. It contains several unique ideas which turn this kit into something special: Two-part, internally cooled trephines allow for a safe and minimal invasive extraction of local bone cylinders. The inner cooling helps to prevent a damage of the bone tissue due to overheating, while the removable trephine working parts allow for an especially easy sampling of the cut out bone cylinders. Either with the aid of the ejection needle or the ejection instruments, which perfectly match with the diameter of the trephine working parts, the bone cylinders can be pushed out of the trephines. Furthermore, the included pre-drilling trephines allow for a safe application of the trephines. Due to the two-part trephines, the variety of different diameters and the optimal matching of the instruments, the Trephine Ejection Kit offers highest efficiency, flexibility and precision to the user for the extraction of bone cylinders.

Art.-No. BTE00

Bone core are harvested from the implant site minimal invasive with the pre- and main trephine in different diameter.

Bone core augmentation stabilized with micro screws.

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
Khoury Karottenheber Khoury Carrot lifter illustrated 1:2	
Ref.	HK001
L ₁ (mm)	38.0
L ₂ (mm)	20.0
L ₃ (mm)	173.0

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AT A GLANCE

- Two-part trephines for a safe and minimal invasive extraction of bone cylinders
- Protection against overheating due to internally cooled trephines
- Innovative application
- Standardized diameters and precise laser markings of the trephines allow for a predictable extraction of defined bone cylinders
- Perfectly matched instruments care for highest precision and flexibility

Micro Screw T5 System Basic NEW

Screw System for Bone Core and Shell Technique |
developed with Prof. Dr. Fouad Khoury

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The Micro Screw T5 Systems are used for the secure fixation and stabilization of bone grafts during bone core or shell technique. The special feature of these systems are the osteosynthesis screws, which are available in two diameters, 1.0 mm and 1.2 mm, and in different lengths. They are made of surgical stainless steel, which gives them a high degree of stability despite their minimized diameter. The new screw drive characterized by a 5-point star-shaped pattern has improved the fit of screwdriver and screw as well as force transmission during insertion.

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Super flexibel | Super flexible

The perfect fit with the new
T5 screw head (torx connection)

	Pilotbohrer Pilot burs		Vorkörner Initial bur		Schraubendreher manuell Screwdriver manual illustrated 1:5		Schraubendreher RA Screw driver RA
Fig.	MSPB1	MSPB2	202RF	202RF	MSTD	MSTDW	
Shank ¹	204	204	206	206	-	204	
Size ²	008	010	008	010	-	-	
Length mm	14.0	14.0	14.0	14.0	148.0	18.0	
mm	0.8	1.0	0.8	1.0	-	-	

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AT A GLANCE

- Safe fixation and stabilization of cortical bone grafts
- Osteosynthesis screws with diameters of 1.0 and 1.2 mm made from surgical stainless steel for high stability
- Self-tapping thread flanks of the screws allow the osteosynthesis screws to be screwed in without creating a gliding hole

	Pilotbohrer Pilot burs		Vorkörner Initial bur		Schraubendreher manuell Screwdriver manual illustrated 1:5		Schraubendreher RA Screw driver RA
Fig.	MSPB1	MSPB2	202RF	202RF	MSTDM	MSTDW	
Shank ¹	204	204	206	206	-	204	
Size ²	008	010	008	010	-	-	
Length mm	14.0	14.0	14.0	14.0	148.0	18.0	
mm	0.8	1.0	0.8	1.0	-	-	

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NEW

Khoury-Line
by
Prof. Dr. Fouad Khoury

Combined System for Extraction and Fixation of Bone Grafts |
developed with Prof. Dr. Fouad Khoury

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
Art.-No. BMTSP

30 Screws

The Micro Screw T5 System Plus is a combined system of the Trephine Ejection Kit, which enables the minimally invasive extraction of local bone cylinders, and the Micro Screw System, which is used for the safe fixation and stabilization of bone grafts. The Micro Screw T5 System Plus thus enables the complete workflow of the Bone-Core-Technique to be carried out. For the extraction of local bone cylinders, the system contains two-part, internally cooled trephines, which, thanks to the removable trephine working parts, enable the bone cylinders to be pushed out particularly easily with the help of special ejection instruments. For the subsequent fixation of the bone cylinders, the Micro Screw T5 System Plus offers osteosynthesis screws with a diameter of 1.0 mm or 1.2 mm, as well as all the necessary instruments for setting the screws. The osteosynthesis screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter. The combination of the two techniques in one system and the optimally matched instruments offer the user a holistic, clear workflow when performing the Bone-Core-Technique and allow for the highest precision from the extraction to the fixation of autologous bone cylinders.

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The perfect fit with the new T5 screw head (torx connection)

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AT A GLANCE

- Combination of systems for carrying out the complete workflow of the Bone-Core-Technique
- Optimally matched instruments for highest precision from extraction to fixation of autologous bone cylinders
- Two-part trephines for the safe and minimally invasive extraction of bone cylinders
- Safe fixation and stabilization of cortical bone grafts
- Osteosynthesis screws with diameters of 1.0 and 1.2 mm made from surgical stainless steel for high stability
- Self-tapping thread flanks of the screws allow the osteosynthesis screws to be screwed in without creating a gliding hole

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
The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

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Super flexibel | Super flexible

The former design of the Khoury
Screw with a square screw head

	Pilotbohrer Pilot burs		Vorkörner Initial bur		Schraubendreher einfach Screwdriver basic		Schraubendreher RA Screwdriver RA	
Fig.	MSPB1	MSPB2	202RF	202RF	MSSDE**		MSSDW**	
Shank ¹	204	204	206	206	-		204	
Size ²	008	010	008	010	-		-	
Length mm	14.0	14.0	14.0	14.0	101.7		24.0	
+ mm	0.8	1.0	0.8	1.0	-		-	

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Super flexibel | Super flexible

The former design of the Khoury Screw with a square screw head

	Pilotbohrer Pilot burs		Vorkörner Initial bur		Schraubendreher manuell Screwdriver manual		Greifer für Schraubendreher manuell Claw for manual screwdriver	Schraubendreher RA Screwdriver RA
Fig.	MSPB1	MSPB2	202RF	202RF	MSSDM**	MSGSD	MSSDW**	
Shank ¹	204	204	206	206	-	-	204	
Size ²	008	010	008	010	-	-	-	
Length mm	14.0	14.0	14.0	14.0	142.2	33.0	24.0	
* mm	0.8	1.0	0.8	1.0	-	-	-	

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
Ost-Tray

Selection of Manual Instruments for Oral Surgery | developed with Prof. Dr. Fouad Khoury

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The Bone Management® Ost-Tray offers a perfect combination of specially configured manual surgical instruments. In addition to a well-thought-out selection of manual instruments, such as high-quality periodontal probes, curettes, scalpel blade holders, mouth mirrors, graduated mouth wedges and various scissors, needle holders, tweezers, clamps, forceps, retractors and tongue depressors, the kit contains an individually designed, strong Khoury-probe as well as raspatories and sharp spoons according to Prof. Khoury.

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


Sinus-Tray

Selection of Manual Instruments for external Sinus Lift | developed with Prof. Dr. Fouad Khoury

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
The new Bone Management® Sinus-Tray is ideal for performing the external sinus lift. The manual instruments designed with Prof. Khoury support the user in gently lifting the Schneiderian membrane and inserting the bone substitute material.

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Inklusive: Bur Block BSCM0 (for Handpiece, 10 slots)

Handinstrumente
Manual instruments


Fig.	MI135	MI136	MI140	MI141	MI142
Name	Kondenser 'Khoury' Condenser 'Khoury'	Knochenmateriallöffel 'Khoury' Bone Material Spoon 'Khoury'	Sinus-Lift Elevator 'Khoury'	Sinus-Lift Elevator 'Khoury'	Sinus-Lift Elevator 'Khoury'
Length mm	168.0	177.0	185.0	190.0	185.0


Soft-Tissue-Tray

Selection of Manual Instruments for Soft-Tissue Management | developed with Prof. Dr. Fouad Khoury

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The new Bone Management® Soft-Tissue-Tray was developed with Prof. Khoury and is perfectly applicable in the context of soft tissue management. A special feature is the scalpel blade holder with swivel head included in the set, with which the scalpel angle can be freely adjusted depending on the processing area. This simplifies handling, especially in areas that are difficult to access.

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
The Smartscraper combines the collection of autologous bone chips and their application to the defect site in a single instrument. The blade produces bone chips, which can then be applied via the syringe mechanism, for example during a sinus lift. The instrument is sterile-packed and ready for immediate use. Maximum bone chip volume: 0.3 cm³ (3 pieces per package, single-use).

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Thanks to its special microblade and slim shape, the Micross bone scraper is particularly suitable for hard-to-reach harvesting sites. Post-operative pain can be minimized as the instrument can be inserted into a tissue tunnel thanks to its small diameter. Maximum bone chip volume: 0.25 cm³ (1 piece per package, single-use).

© Meda

Smartscraper
Smartscraper

Fig.

BS005

Micross
Micross

Fig.

MBHD1


SAFESCRAPER® Twist

Ready-To-Use Cortical Bone Collector

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With the Safescraper® Twist, autologous cortical bone chips can be obtained easily and safely. These are safely collected in the transparent chamber while the Safescraper® is in use - for optimal control of the amount of bone harvested. The bone material obtained can be used immediately thanks to the immediate mixture with blood and is suitable for filling up defects of all kinds.

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AT A GLANCE

- Ready to use
- Excellent cutting performance for a minimally invasive and patient-friendly procedure
- Straight and angled design for every situation
- Half-round blade with 160° cutting area for maximum efficiency on all bone surfaces
- Capacity of the chambers of 3.0 cm³ (straight version) or 2.5 cm³ (angled version)
- Safety mechanism prevents the chambers from opening during bone harvesting

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The Tenting Screw System contains a broad selection of screws and instruments for use in intraoral bone augmentation and offers an alternative to the shell technique and the use of titanium-reinforced PTFE membranes.

The tenting screws are used for the fixation and stabilization of bone grafts and bone substitutes as well as for the mechanical support of barrier membranes. The screws contained in the system are made of the titanium alloy Grade 5-ELI (Titanium Grade 23) and have a self-tapping thread. The screw head includes an interface for holding the insertion tool and for optimal force transmission (hexalobular socket).

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AT A GLANCE

- System for the treatment of simple and complex bone defects by means of screw-assisted regeneration techniques
- Alternative to shell technique and the application of titanium-reinforced PTFE membranes
- Reliable fixation and stabilization of bone blocks as well as bone substitute materials
- Efficient insertion of tenting screws due to self-cutting thread design
- High biocompatibility due to titanium alloy

	Vorkörner Initial bur	Spiralbohrer Twist drills	Schraubendreher RA lang Screwdriver RA long
Fig.	186RF	203RF	31053
Shank ¹	204	205	205
Size ²	018	011	-
Length mm	12.0	9.0	27.0
* mm	1.8	1.1	-

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The Umbrella Screw is intended for use in intraoral bone augmentation and offers an alternative to the shell technique and the use of titanium-reinforced PTFE membranes. The Umbrella Screw is used for the fixation and stabilization of bone grafts and bone substitutes as well as for the mechanical support of barrier membranes. The screw is made of the titanium alloy Grade 5-ELI (Titanium Grade 23) and has a selftapping thread. The screw head includes an interface for holding the insertion tool and for optimal force transmission (hexalobular socket).

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AT A GLANCE

- Alternative to shell technique and the application of titanium-reinforced PTFE membranes
- Reliable fixation and stabilization of bone blocks as well as bone substitute materials
- Efficient insertion of tenting screws due to self-cutting thread design
- High biocompatibility due to titanium alloy

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MEMBER OF THE GLOBAL
BONE MANAGEMENT® &
RESTORATIVE NETWORK

Developed with
PD Dr. Dr. Florian Bauer
Miesbach, Germany

Praxisklinik MKG Miesbach
www.mkg-miesbach.de

My First GBR Kit NEW


Selection of Manual Instruments for GBR

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The kit consists of 12 manual instruments selected from our portfolio for application in guided bone regeneration techniques (GBR). It is optimally suited for young surgeons who are taking their first steps in regenerative surgery. The instruments can be safely transported and sterilized within the tray.

Art.-No. GBR00

Handinstrumente
Manual instruments


Fig.	MI070	MI071	MI077	MI139	MI057	MI101
Name	Raspatorium Molt Periosteal Elevator Molt	Raspatorium Khoury Periosteal Elevator Khoury	Skapellklingenhalter Scapel Handle	Kürette Columbia Curette Columbia	Scharfer Löffel 3 mm Bone Curette 3 mm	Wundhaken Langenbeck 42x10 mm Tissue Retractor Langenbeck 42x10 mm
Length mm	180.0	183.0	152.0	166.0	175.0	220.0

Handinstrumente
Manual instruments

Fig.	MI089	MI106	MI041	MI048	MI050	MI072
Name	Chirurgische Schere Kelly Surgical Scissors Kelly	Hohlmeißelzange Blumenthal Bone Rongeur Forceps Blumenthal	Nadelhalter Crile Wood, TC Needle Holder Crile Wood, TC	Pinzette Tweezers	Pinzette, TC Tweezers, TC	Raspatorium U Periosteal Elevator U
Length mm	160.0	155.0	150.0	175.0	175.0	172.0

Transfer-Control

Horizontal and Vertical Bone Replacing System

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Transfer-Control permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in five diameters. They are perfectly adapted to each other so that the outer diameter of the ablative bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm and 8.0 mm.

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

AT A GLANCE

- Easy and safe transplantation of bone cylinders
- Intelligently coordinated instruments
- Perfectly fitting bone cylinders
- Precise work due to depth marks
- Optimally suitable for vertical and horizontal augmentations
- Creation of refreshed contact surfaces for rapid vitalization and wound healing

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The most important condition for a safe healing of a bone graft is a congruent and fresh contact surface of the implant area. The Transfer-Ring-Control I System allows for this precondition for the healing of vertical bone grafts in a simple and controlled manner. This system offers a variety of various instrument sizes and provides high flexibility for the extraction of bone rings. Depending on requirements, various sizes of Initial Bur Trephines, Trephines and Ablative Burs are available. With the aid of the additional tungsten carbide bur, sharp bone edges can be smoothed precisely.

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Transfer-Ring-Control II

One-Stage Bone Ring Augmentation System | Technology by Dr. Bernd Giesenhagen

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Transfer-Ring-Control II is a system for the bone ring technique, developed by Dr. Giesenhagen. It allows for one-stage proceeding of augmentation and implantation in case of resorbed jaw bone, if implant insertion alone is not possible because of inadequate spatial conditions. Due to the perfectly matched instruments, there is created a congruent bone and implant site, which allows for the treatment of three-dimensional defects with the help of prefabricated, allogenic bone rings. Because of the optimal fit between jaw bone, bone ring and inserted implant and the possibility of the one-stage proceeding, high primary stability is accomplished and healing time is significantly reduced.

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


Technology by
Dr. Bernd Giesenhagen
Kassel, Germany

Transfer-Control Plus

Combined System for Horizontal and Vertical Bone Augmentation

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
Transfer-Control Plus permits precise and standardized transplantation of bone cylinders and rings for horizontal and vertical bone augmentation. Different types of bone burs allow for high flexibility for the extraction of bone grafts. The instruments are perfectly adapted to each other often producing a press fit between the cylindrical bone or ring and the prepared recipient site. This press fit can then be strengthened with a fixation screw which is sold separately (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound-healing. So a transplantable bone site is reached after only 3-4 months.

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AT A GLANCE

- Combined system for easy and safe transplantation of bone cylinders and rings
- Intelligently coordinated instruments
- Optimally suitable for vertical and horizontal augmentations
- Precise work due to depth marks
- Creation of refreshed contact surfaces for rapid vitalization and wound healing

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

The Osteosynthesis System Screw System TX enables the universal use of fixation for autologous bone cylinders, bone blocks and bone shells for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

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
The Osteosynthesis System Screw System TX Professional enables the universal use of fixation for autologous bone cylinders, bone blocks and bone shells for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

The Screw System TX Professional contains additional screws in the lengths of 13.0 mm and 16.0 mm.

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

The Trephine Basic Kit provides a selected set of high-quality trephine drills in different sizes. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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AT A GLANCE

- Selection of high quality trephines for the extraction of bone cylinders
- Various trephine sizes for every application
- Precise depth markings for optimal control
- High cutting performance for good handling

Trephine Kit

Trephine System for Extraction of Cylindrical Bone Grafts

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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AT A GLANCE

- Broad selection of high quality trephines for the extraction of bone cylinders
- Various trephine sizes for every application
- Precise depth markings for optimal control
- High cutting performance for good handling

Trephine Kit Long

Trephine System for Extraction of Cylindrical Bone Grafts

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 14 mm* and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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Punch Basic Kit

Biopsy Punch Set for Extraction of Tissue Grafts

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The Punch Basic Kit contains Biopsy Punches in five different diameters for the performance of tissue removal. They serve for the application in immediate implant placement, for taking mucosa grafts in the scope of soft tissue management or biopsy as well as for uncovering the gingiva after implantation. The Biopsy Punches create clean cut surfaces and care for minimal tissue traumatization.

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Saw Basic Kit

Selection of Various Diamond Coated Saws for Oral Surgery

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The Saw Basic Kit provides a selection of diamond saws in the three most important diameters. They are used for the application in bone spreading or splitting operations. Due to the various sizes, the set offers the clinician the possibility of having a proper selection of these saws during surgery.

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© Dr. Beyer & Dr. Kistler

Crestal-Lift-Control Basic

Crestal Approach Sinus Lift System

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Crestal-Lift-Control Basic is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control Basic is an introductory system with reduced instrument set.

CAUTION: To prevent injury of any anatomical structure, the apical extra lengths of the CL001 of 0.58 mm must be considered.

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
Crestal-Lift-Control is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control provides an extensive range of instruments for the performance of a safe internal sinus lift.





CAUTION: To prevent injury of any anatomical structure, the apical extra lengths of the CL001 of 0.58 mm must be considered.

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External-Lift-Control is a system for the easy performance of an external sinus lift. After lateral access through the facial wall of the maxillary sinus, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. Then, a suitable augmentation material can be inserted under the membrane. The sandblasted silver-matt hand instruments provide a reflection-free working environment.

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
Internal-Lift-Control is a system for the performance of an internal sinus lift. With the aid of the elevators, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 5 mm.

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Lift-Control Plus

Combined System for Internal and External Sinus Lift

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The combined Lift-Control Plus system provides perfectly matched instruments for the performance of internal as well as external sinus floor elevation. For the internal sinus lift, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 5 mm. After lateral access through the facial wall of the maxillary sinus for the performance of the external sinus lift, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. After this, a suitable augmentation material can be inserted under the membrane.

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Surgical Kit 1

Lateral Approach System for External Sinus Elevation

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The Surgical kit 1 contains instruments with handpiece shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. The large diameter of the instruments protects the membrane from unintentional perforation.

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AT A GLANCE

- Compilation of handpiece instruments for the lateral approach during the external sinus lift
- Round instruments protect the Schneiderian membrane from unintentional perforation
- Ideal for gradual, safe thinning of the facial wall of the maxillary sinus
- Round diamonds and carbide burs in three different diameters

Implant Basic Kit NEW

Selection of Manuel Instruments for Implantology |
developed with Prof. Dr. Georg-Hubertus Nentwig

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The Bone Managment® Implant Basic Kit offers a perfect combination of specially configured manuel surgical instruments for implantology. The kit contains two different raspatories, a scalpel handle, surgical tweezers, two different probes, a bone curette, a haemostatic forceps, a needle holder, curved scissors and a cheek retractor.

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AT A GLANCE

- Designed for surgical interventions in the context of implantology
- Surgical instruments can be prepared and transported within the tray
- Instrument composition according to Prof Hubertus Nentwig, the developer of the myplant implant

Developed with
**Prof. Dr. Georg-Hubertus
Nentwig**
Germany

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The Surgical Kit 3 contains conical Tungsten Carbide Burs with friction grip shank of various sizes and toothing which allow for a broad field of application. The particularly easy cutting and efficient cutters have an extra-long shank for an optimal view on the working area and sharp tips for the safe application in bone.

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The Surgical Kit 4 contains a selection of conical Tungsten Carbide Cutters with handpiece shank for a broad field of application in oral surgery. The cutters have a particularly easy cutting and efficient saw toothing and a specifically designed twist angle. Further, the narrow working parts allow for an optimal view on the working area.

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* 

Alveoplasty Surgical Kit

Surgical Contouring and Modification Kit

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
Art.-No. ALV18

The Alveoplasty Surgical Kit contains the most important instruments for surgical contouring and modification of the alveolar ridge when unwanted unevenness of the jawbone occurs after tooth extraction or bone degeneration. In this case, the kit offers the best way to smooth the jawbone. The Alveoplasty Kit contains two initial burs and two Lindemann burs, which allow for the performance of precise, initial osteotomies, even in hard cortical bone. To perform the alveoplasty itself, five different carbide cutters and a round diamond are included. The round working parts optimally protect the soft tissue against injuries. Thus, the alveolar ridge is prepared carefully and precisely for implant insertion or denture restoration.

	Vorkörner (2 Stück) Initial burs (2 pcs.)	Chirurgiefräser Surgical Cutters		Hartmetallinstrumente Tungsten Carbide Instruments					Diamantinstrument Diamond Instrument	Kieferkammkonturierer Ridge Contouring Bur
Fig.	186RF	HM166**	165RF	HM71**	HM71**	HM72GX**	HM77GX**	HM251GX**	801*	RCB00
Shank ¹	204	205	205	104	104	104	104	104	204	205
Size ²	018	021	023	040	050	060	060	040	033	075
Length mm	12.0	11.0	7.0	-	-	12.0	11.5	8.3	-	6.0
Δ mm	-	-	-	-	-	-	-	-	-	4.2
+ mm	1.8	2.1	2.3	4.0	5.0	6.0	6.0	4.0	3.3	7.5

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Guided-Drill-Stop-Control

Guided Drill System for Safe Implant Site Preparation | developed with PD Dr. Dr. Florian Bauer

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

The Guided-Drill-Stop system is used for the easy and safe preparation of the implant site. The special feature of the system are the specially developed twist drills with guiding tip, with which an intermediate step is added prior to the actual extension drilling. With this step, the upper part of the implant site is expanded to the desired diameter, while the lower part is used to guide the drill. This avoids deviation from the desired axis. In addition, the stop sleeves provide the defined depth stop during all drillings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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
Screw and sleeve for safe intraoral fixation of surgical guides, screw self-drilling and self-tapping, manual insertion or by means of screw driver (e.g. TL0T1).

GBMN
MEMBER OF THE GLOBAL
BONE MANAGEMENT® &
RESTORATIVE NETWORK

Developed with
PD Dr. Dr. Florian Bauer
Miesbach, Germany

Praxisklinik MKG Miesbach
www.mkg-miesbach.de

The K Kit Implantoplasty NEW

Tungsten Carbide Finishers for Periimplantitis Treatment | developed with Dr. Philippe Khayat

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The kit comprises five carbide finishers and a high-quality diamond instrument. All instruments have been assembled in one kit to enable precise and effective smoothing of rough and contaminated implant surfaces. By smoothing the exposed part of the implant and treating the area appropriately, it enables implant preservation and therefore minimally invasive treatment of periimplantitis. The extra-long shank on all instruments ensures easy handling in hard-to-reach areas, while the different working part shapes provide an optimized equipment for smoothing of all types of implant geometries.

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


AT A GLANCE

- Matching carbide finishers and one diamond bur for effective intraoral smoothing and cleaning of contaminated implant surfaces by means of implantoplasty
- XXL FG shank facilitates work on hard-to-reach implants
- Different shapes and sizes of instruments for each implant and patient situation

MEMBER OF THE GLOBAL
BONE MANAGEMENT® &
RESTORATIVE NETWORK

Developed with
Dr. Philippe Khayat
Paris, France

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The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the idea peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothing. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access, deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counter-clockwise so the implant surface is smoothed.

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AT A GLANCE

- Perfectly matched tungsten carbide finishing burs for the effective, intraoral treatment of titanium
- XXL FG shanks facilitate the treatment of difficult to access implants
- Different shapes and sizes for every individual patient case
- Clean and smooth titanium surfaces

Degranulation Kit

Diamond System for effective Bone Debridement

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The Degranulation Kit contains four extra-coarse, spherical diamond instruments in diameters of 1.0, 2.5, 3.0 and 3.5 mm. In general, this kit is suitable for removing unwanted residues in extraction sockets. In particular, however, the diamond instruments have been developed to enable quick and easy bone debridement. With the help of the diamonds, degranulation tissue adhering to the bone can be removed quickly and easily without damaging the surrounding bone, so that subsequently inserted augmentation material comes into direct contact with healthy bone. The grooves on the shaft of the diamond instruments prevent blood from contaminating the contra-angle handpiece holder.

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AT A GLANCE

- Extra coarse, round diamonds for quick and easy debridement
- Reliable removal of degranulation tissue without damaging the surrounding bone
- Four different diameters for use in all situations
- Extra long instrument necks for an optimal view of the operating area
- Optimal addition to the degranulation process


MEISINGER

dPTFE Membrane

Micro-textured, high-density PTFE Membrane

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MEISINGER's dPTFE membranes are made of high-density PTFE with a pore size of less than $0,3 \mu\text{m}$. Thus, the non-resorbable membrane is insensitive to bacteria and allows a maximum of predictability and aesthetics. When using the dPTFE membranes, no primary wound closure is necessary. The exposed membrane can be removed with tweezers after at least 21 days without surgical intervention, no additional anesthesia is required. The hexagonal depressions are designed to increase the surface area, which leads to an increased stability of the membrane and favors cell adhesion without increasing the porosity. For best results, the side with the depressions should be placed in the direction of the soft tissue. The Meisinger dPTFE membranes can be cut to size as needed.

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AT A GLANCE

- Non-resorbable membrane made of high-density PTFE
- Maximum of predictability and aesthetics
- No primary wound closure necessary
- Non-surgical removal of the membrane
- Hexagonal depressions of the surface for high stability
- No disruption of wound healing due to early resorption

MEISINGER Titanium-Reinforced dPTFE Membrane

Titanium-reinforced, high-density PTFE membrane

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MEISINGER's Titanium-Reinforced membranes are made of high density PTFE. In combination with titanium reinforcement, the membranes offer exceptional dimensional stability and optimum flexibility at the same time. Therefore, the membranes are particularly suitable for the treatment of 3- to 4-wall defects in the course of the alveolar ridge augmentation. The wide variety of membrane shapes and sizes covers a variety of different defects and minimizes the need for trimming. If necessary, however, the membranes can be optimally trimmed to the required shape. The titanium-reinforced membranes are available in thicknesses of 150 μm and 250 μm .

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AT A GLANCE

- Non-resorbable membrane made of high-density PTFE
- Fine Titanium Reinforcement from Titanium Grade I
- Optimum combination of flexibility and dimensional stability
- Treatment of 3- to 4-wall defects
- Wide range of different membrane shapes
- Availability of two different membrane thicknesses

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Anwendung bei großen bukkalen Defekten.
Designed for large buccal defects.

Art.-No. MTi2502536-1 Art.-No. MTi2502536-2
Art.-No. MTi1502536-1 Art.-No. MTi1502536-2
25 mm x 36 mm 25 mm x 36 mm
(1 membrane / box) (2 membranes / box)

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
The unique circular macroporous design of the Reinforced PTFE Meshes allows for direct contact between the bone graft and periosteum, allowing naturally occurring revascularization and infiltration of cells into the bone graft. The fine titanium frame maintains space essential for horizontal and vertical ridge augmentation. The PTFE Meshes have a thickness of 200 µm and easily conform to tissue contours.

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AT A GLANCE

- Non-resorbable PTFE Mesh
- Circular macropores for direct contact between bone graft and periosteum
- Fine titanium frame maintains space for horizontal and vertical ridge augmentation
- Maximum of predictability and aesthetics
- No disruption of wound healing due to early resorption

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Art.-No. RPM200PDML
(Unterkiefer links / Mandible left)
38 mm x 38 mm
(1 membrane / box)

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
MEISINGER's bovine collagen membrane is made from highly purified type 1 collagen derived from bovine achilles tendon. A great advantage of the bovine membranes is their ideal handling. They combine high flexibility with optimum stability. This makes them flexible, adaptable to any shape and easy to place. At the same time they are stable enough to hold their shape and not to collapse into the defect. The unique fiber orientation ensures maximum tensile strength. Thus, the membranes can be safely fixed by sutures or pins without tearing. Due to the long and predictable resorption time of 26 to 38 weeks, there is no risk of particle loss due to early resorption.

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AT A GLANCE

- Long resorption time of 26 to 38 weeks
- Ideal handling
- High flexibility for easy ridge preservation
- Optimum dimensional stability
- High tensile strength for safe fixing

myAligner

Meisinger

THE ALIGNER SYSTEM

MADE FOR DENTISTS AND
PRE-SURGICAL ALIGNER
THERAPIES

NEW!

CLASSIC

Biomechanically
optimized thermo-
formed aligners

3D

Direct printed
shape memory
aligners

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Zmatrix™ is a natural, native collagen membrane. Cross-linking chemicals and agents are unnecessary. Proprietary processing technology allows preservation of collagen as well as extracellular components including laminin, fibronectin, elastin, and glycosaminoglycans.*

Art.-No. ZM1520
15 mm x 20 mm
(1 membrane / box)

Art.-No. ZM2530
25 mm x 30 mm
(1 membrane / box)

Art.-No. ZM3040
30 mm x 40 mm
(1 membrane / box)

*Hoganson DM, Owens GE, O'Doherty EM, Bowley CM, Goldman SM, Harilal DO, Neville CM, Kronengold RT, Vacanti JP. Preserved extracellular matrix components and retained biological activity in decellularized porcine mesothelium. Biomaterials. 2010; 27: 6934-6940.

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AT A GLANCE

- Perfectly soft consistency
- Easy to handle - Designed to drape without adhering to itself
- Elastic - Natural peritoneum collagen structure allows for elasticity
- Resorption time of about 26 weeks

Cytoplast™ RTMPlug, RTMFoam, RTMTape

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Note: Wound dressings will be essentially resorbed within 30 days.

RTMPLUG

Art.-No. RTMPLUG10
10 mm x 20 mm
(10 / box)

RTMFOAM

Art.-No. RTMFOAM10
20 mm x 40 mm
(Dicke/Thickness: 3 mm, 10 / box)

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APPLICATIONS:

- Surgical wounds
- Periodontal surgical wounds
- Extraction sites
- Dental sores
- Oral ulcers (non-infected or viral)
- Suture sites
- Burns
- Traumatic wounds

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

The high quality of the MEISINGER PTFE suture allows for atraumatic wound closure after dental procedures. The high-quality needles are made of stainless steel and have a very high bending strength with a minimum of wear. Sharp needles are the result, even after repeated piercing. The constant cutting performance allows for a particularly tissue-saving work. The suture is made of 100% medical grade bioinert PTFE. Thanks to its smooth surface, the monofilament prevents the adhesion of bacteria to the thread and ensures a resistance-free gliding through the tissue. The PTFE is a soft, non-resorbable material with consistent tensile strength, allowing the user to work very precisely, planning with a defined healing time, and patient-friendly care.

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

Suture: USP 3-0

Needle: 16 mm, 3/8 Circle, reverse cutting, black coated

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Suture: USP 3-0

Needle: 19 mm, 3/8 Circle, reverse cutting, black coated

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
The Pro-Fix™ Precision Fixation System is manufactured to precise tolerances to ensure easy pick-up of screws, stable transfer to the surgical site and quick engagement in cortical bone. All kit components are labeled, organized, stored and sterilized together to simplify things for the surgeon and surgical team.

Tecapro™ Storage Tray, Tecapro™ Storage Tray

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The autoclavable Tecapro™ Storage Tray is cleverly designed to store not only fixation instruments, but also a variety of membrane fixation, bone fixation, and tenting screws needed in dental bone grafting.

Schrauben-Organizer, Screw Organizer Dial

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Up to 100 fixation screws – including bone fixation, membrane fixation, and tenting – may be stored in a built-in screw organizer dial designed and labeled for easy identification, easy storage and simple re-ordering.

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The stainless steel driver handle and cruciform driver blade fit together to make up a drive system that ensures easy pick-up of screws and stable transfer to the surgical site.

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Blades work universally with all Pro-Fix™ screws.

Art-No.	Inhalt/Content	
PFT	Autoklavierbares Tecapro™ Storage Tray mit Schrauben-Organizer Autoclavable Tecapro™ Storage Tray with screw organizer dial	
PFDH	Schraubendreher-Griff aus rostfreiem Stahl Stainless steel driver handle	
PFDB	76 mm Kreuzschraubendreher 76 mm cruciform driver blade	
PFDB56	56 mm Kreuzschraubendreher 56 mm cruciform driver blade	
BI1001	1,2 mm Durchmesser RA-Pilotbohrer 1.2 mm diameter latch type pilot drill	
PFDBCA	RA-Schraubendreher (24 mm lang, 10 mm exponierte distale Länge) Contra angle blade (24 mm long, 10 mm exposed distal length)	

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Pro-Fix™ Membrane Fixation Screws are designed as an attractive alternative to using tacks for membrane stabilization. Easy pick-up, solid stability of the screw during transfer to the surgical site, and easy placement make membrane fixation fast and easy.

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Further self-drilling Membrane Fixation Screws

Art-No.	verpackt à packed à	Inhalt/Content	
PFMF-5*	5	1,5 mm x 3,0 mm selbstschneidende Schrauben aus Titanlegierung zur Membranfixierung 1.5 mm x 3.0 mm self-drilling titanium alloy membrane fixation screws	

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Pro-Fix™ Tenting Screws are designed with a self-drilling tip, polished neck and broader head to maintain space under resorbable and non-resorbable membranes in horizontal and vertical bone regeneration procedures.

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Further self-drilling Tenting Screws

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Pro-Fix™ Bone Fixation Screws are designed with finer pitched, self-tapping threads that give the screws greater clamping force while using less driver torque. The screws' threads are equipped with a cutting flute that allows for easier insertion into harder bone. The screws are placed into a 1.2 mm pre-drilled pilot hole.

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Further self-tapping Bone Fixation Screws

Art-No.	verpackt à packed à	Inhalt Content	
PFB8*	1	1,5 mm x 8,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 8.0 mm self-tapping Bone Fixation Screw	
PFB8-5*	5	1,5 mm x 8,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 8.0 mm self-tapping Bone Fixation Screws	
PFB10*	1	1,5 mm x 10,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 10.0 mm self-tapping Bone Fixation Screw	
PFB10-5*	5	1,5 mm x 10,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 10.0 mm self-tapping Bone Fixation Screws	
PFB12*	1	1,5 mm x 12,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 12.0 mm self-tapping Bone Fixation Screw	
PFB12-5*	5	1,5 mm x 12,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 12.0 mm self-tapping Bone Fixation Screws	
PFB14*	1	1,5 mm x 14,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 14.0 mm self-tapping Bone Fixation Screw	
PFB14-5*	5	1,5 mm x 14,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 14.0 mm self-tapping Bone Fixation Screws	

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
cerabone® plus is a combination of cerabone® granules and sodium hyaluronate, which is a water-soluble salt form of hyaluronic acid. Thanks to the pronounced liquid binding capacities of hyaluronate, cerabone® plus, upon hydration, forms a malleable material that enables both easy uptake and delivery to the site of application. cerabone® plus provides a handling alternative by combining the established bone grafting material cerabone® with the well-known properties of hyaluronic acid.

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AT A GLANCE

- Sticky and malleable following hydration
- Efficient defect filling and time-saving application
- Easy defect contouring
- Minimized displacement of single granules during application



collafleece®

Resorbable collagen for wound management

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collafleece® is a wet-stable, porcine collagen sponge with highly efficient hemostatic properties. Its natural porous collagen structure supports the hemostasis and controls the natural healing of the wound. Owing to its spongy collagen structure, the collafleece® resorbs within 2-4 weeks.

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AT A GLANCE

- Highly efficient hemostatic collagen
- Fast resorption by enzymatic degradation (2-4 weeks)
- Easy application
- Maintains integrity in presence of blood and during application
- Wound protection / wound healing

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NanoBone® Catalogue
separately available

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AT A GLANCE

Synthetic bone grafting material

NanoBone® contains absolutely no human or animal ingredients so no material-related contamination risks and alleviates ethical concerns for patient consent. High-tech production processes ensure consistent and excellent quality.

Complete remodelling

Thanks to its special structure, NanoBone® can be used for natural remodelling. Osteoclasts resorb the material as osteoblasts form new bone. Material resorption and bone formation are linked. This means predictable results and a reliable bone basis.

Controlled osteoinduction

NanoBone® is osteoconductive and osteoinductive. In other words, it is not simply a scaffold: it actively promotes bone formation.

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NanoBone® | granulate is available in either fine or coarse granules to suit the indication.

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The granulate can be mixed with the patient's blood or with a sterile saline solution. The manufacturer recommends mixing with blood. The material is extremely hydrophilic and fully absorbs the blood or saline solution quickly.

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When mixed with blood, **NanoBone® | granulate** takes on a paste-like consistency and can be easily applied with a spatula or augmentation spoon.

NanoBone® | QD

Ideal bone healing in combination with perfect handling

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AT A GLANCE

The applicator


- "Ready to use" direct application from the applicator
- Easy, intuitive use
- Precise placement of the **NanoBone® QD** material
- Controlled and exact portioning
- Allows for small approaches to the application site (outer diameter of the applicator just 8 mm)
- Allows for retrograde filling of bone defects
- Difficult to access bone defects treated due to the long applicator

The material

- Shapeable and adaptable to the defect
- Strong cohesion of the material
- Sticks to the bone (adhesion)
- Stable also in aqueous environment and profusely bleeding wounds
- Controlled osteoinduction

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Easy-Clean Scissors

Surgical Scissors for easy and safe processing

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The MEISINGER Easy-Clean Scissors made of high-quality stainless steel can be easily and completely disassembled and put back to-gether again. Residues below the hinge can be completely removed – for a residue-free processing*. For dismantling, the scissors only need to be opened at a right angle. Then the upper part of the scissors is lifted briefly. In addition, the scissors have a black coating that is applied using PVD procedure (vacuum-based coating). Thanks to this black coating the products reach high hardness and durability. Furthermore, disturbing light reflections are reduced, and glare-free work is possible.

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AT A GLANCE

- Surgical scissors made of high-quality stainless steel
- Easy and complete dismantling for residue-free processing*
- Black coating for high hardness and durability as well as glare-free work
- Scissors in two different versions

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MEISINGER offers one of the most extensive programs in the area of the Surgical Round Drills, Surgical Cutters, and Trephines. The optimized special construction of the instruments guarantee a safe application and fast achievement for the desired working result.

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Surgical Cutters of Tungsten Carbide

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The special, easy-cutting and efficient saw toothing and the specifically selected angle of twist allow a wide field of application of the instruments within the field of oral surgery.

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Surgical round drill of tungsten carbide

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
Surgical round drill of stainless Steel

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
Ridge Contouring Bur of stainless steel

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

* im Saw Basic Kit 7150 enthalten * contained in the Saw Basic Kit 7150

915DM mittel / medium

Fig.	Schaft / Shank	Ref.-Nr. / Ref.-No.		1
			L mm	0,25
915DM	HP	806 104 35M 524		220
	unmounted	806 900 35M 524		220

Super flexibel | Super flexible

Vorkörner

Initial Burs

186RF Edelstahl / stainless steel

Fig.	Schaft / Shank	Ref.-Nr. / Ref.-No.		5
			L mm	12,0
186RF	HP	330 104 684 377		018
	RA	330 204 684 377		018

Vorkörner, Stahl rostfrei | Initial bur, stainless steel

187RF Edelstahl / stainless steel

Fig.	Schaft / Shank	Ref.-Nr. / Ref.-No.		5
			L mm	19,0
187RF	RA	330 204 685 377		018

Vorkörner mit Stop, Stahl rostfrei | Initial bur with Stop, stainless steel

188RF Edelstahl / stainless steel


Fig.	Schaft / Shank	Ref.-Nr. / Ref.-No.		2	2
			L mm	17,3	27,3
188RF	RAXL	330 206 686 378		014	
	RAXXL	330 207 686 378			014

! Developed with Prof. Dr. Fouad Khoury

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



Innengekühlter Trepan-Schaft | Internally cooled trepan shank

229KH*


Schaft / Shank	Fig.	229KH	229KH	229KH	229KH
unmontiert / unmounted		021	025	029	033
L mm		15,5	15,5	15,5	15,5
1		3,1	3,5	3,9	4,3
2		2,1	2,5	2,9	3,3
3		1	1	1	1

Innengekühlte Trepan-Arbeitssteile | Internally cooled trepan work parts

! Technology by Dr. Bernd Giesenhausen

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Stainless Steel Bur Blocks, non-corrosive, sterilizable

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MEISINGER drill stands are ideally suited for the effective cleaning, sterilization and organization of the required instruments.

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The following reference values for rotation speeds apply to surgery in general:

Hand piece (HP):

Optimum: 6.000 - 10.000 rpm
Maximum: 40.000 - 50.000 rpm

Angle piece (RA):

Optimum: 6.000 - 10.000 rpm
Maximum: 40.000 - 50.000 rpm

Friction Grip (FG):

Optimum: 80.000 rpm
Maximum: 100.000 - 120.000 rpm

Please also observe adapting the rotation speed in relation to the diameter of the instrument as well as the prevailing indication and enough cooling. As a general rule, the larger the working part of an instrument, the lower the speed should be set.

Please also pay attention to the special preparation specifications for products made of tool steel.

Please also note the recommended and maximum rotary speeds for the individual instruments on the product packaging and in the corresponding instructions for use.

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The symbols give merely suggestions for the possible implementation of the products. The user decides and takes full responsibility about the precise deployment according to existing indications. Please follow general application and safety instructions for MEISINGER products in the medical and dental area and also the advice for processing. Details can be found on the internet under en.meisinger.de/services/downloads/ or you can request one by mail.

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Cleaning/Disinfection/Sterilization

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Contains hazardous substances: Cobalt, CAS: 7440-48-4

This product contains cobalt in more than 0.1 percent by mass and is therefore, in accordance with current regulations, to be labeled as CMR material class 1B (carcinogenic, mutagenic and/or reprotoxic). It has been proven that there is no increased risk of cancer or adverse reproductive or genetic effects when the product is used according to the directions for its intended use.

Art. Nr.	Name	Page
2575	Periimplantitis Kit	90
2580	The K Kit Implantoplasty	88-89
7120	Trephine Basic Kit	70
7121	Trephine Kit	71
7122	Trephine Kit Long	72
7140	Punch Basic Kit	73
7150	Saw Basic Kit	73
31054	Schraubendreher RA kurz / Screwdriver RA short	68
203S	Spiralbohrer mit Stop zum Dekortikalisieren / Twist Drill with Stop for decorticating	26, 27
915DM	Diamantierte Scheibe / Diamond Disc	44, 48, 49, 54
943DC	Diamantierte Scheibe / Diamond Disc	54
ALV18	Alveoplasty Surgical Kit	85
BBR00	Transfer-Ring-Control II	65
BCL00	Crestal-Lift-Control	75
BGS01	Guided-Drill-Stop-Control	86-87
BIB00	Implant Basic Kit	82-83
BIE00	Implant Extraction Kit	41
BIR00	Implant Removal Kit	40
BKM00	Master-Mill	24-25
BLIPL	Lift-Control Plus	78-80
BMCBA	Master-Core Basic	22
BMCPR	Master-Core Professional	23
BMP00	Master-Pin-Control	27
BMPBA	Master-Pin-Control Basic	26
BMPPL	Master-Pin-Control Plus	28-29
BMPXB	Master-Pin-X-Control Basic	30
BMPXU	Master-Pin-X-Control Upgrade	31
BMS00	Micro Screw System	49
BMSBA	Micro Screw System Basic	48
BMST0	Micro Screw T5 System	45
BMSTB	Micro Screw T5 System Basic	44
BMSTP	Soft Tissue Punch System	38
BMTSP	Micro Screw T5 System Plus	46-47
BOST0	Ost-Tray	50-51
BR001	Knochenringpinzette / Bone Ring Tweezers	65
BRSS1	Root Sectioning Kit Samarani 1	12
BRSS2	Root Sectioning Kit Samarani 2	13
BS003	Safescraper® TWIST straight Version	56
BS004	Safescraper® TWIST Curve angled Version	56
BS005	Smartscraper	55
BSIN0	Sinus-Tray	52
BSK02	Surgical Kit 2	81
BSK03	Surgical Kit 3	84
BSK04	Surgical Kit 4	84
BSKSL	Surgical Kit 1	80
BSOT0	Soft-Tissue-Tray	53
BTE00	Trephine Ejection Kit	42-43

Art. Nr.	Name	Page
BTL00	External-Lift-Control	76
BTS00	Tenting Screw System	58
BTX00	Screw System TX	68
BTXPR	Screw System TX Professional	69
CBE00	Benex®-Control	08-09
CCLBA	Crestal-Lift-Control Basic	74
CCR00	Crest-Control	14-15
CL021	Tiefenmess-Sonde / Depth gauge	74
CM001	Easy-Clean Scissors	112-113
CM002	Easy-Clean Scissors	112-113
CSL00	Internal-Lift-Control	77
CSP12	Split-Control 12 mm	16-17
CSPPL	Split-Control Plus	18-19
CTR00	Transfer-Control	62-63
CTR10	Transfer-Ring-Control I	64
CTRPL	Transfer-Control Plus	66-67
DBH06	Führungshülse Drill sleeve	87
DEG00	Degranulation Kit	91
GBR00	My First GBR Kit	60-61
GF080	Surgical Guide Fixation Screw	87
GTAOK	Get 'em all out Kit	10
HF020	Master Surgery Kit	32-33
HK001	Khoury Karottenheber / Khoury Carot lifter	42
KD01G	Diamantinstrument zur Modifizierung des Kieferkamms / Diamond instrument for modifications of the alveolar ridge	85
KD02G	Diamantinstrument zur Modifizierung des Kieferkamms / Diamond instrument for modifications of the alveolar ridge	85
MBHD1	Micros	55
MI154	Hammer / Hammer	77, 79
MI159	Knochenblockpinzette / Bone Block Tweezers	54
MI164	Spülkanüle mit Olive / Irrigation Needle with Bulb End	50, 52
PFBK12	Pro-Fix™ Bone Fixation Kit	105
PFMK20	Pro-Fix™ Membrane Fixation Kit	103
PFTK12	Pro-Fix™ Tenting Kit	104
RTR00	RTRO Kit	11
SL1	Einschlagwerkzeug / Impact tool	79
TL0T1	Schraubendreher manuell / Screwdriver manual	68, 69
USL600	Umbrella Screw	59
/	Membranen / Membranes	92-100
/	Nahtmaterial / Suture	101
/	cerabone® plus	106
/	collafleece®	107
/	NanoBone®	108-111
/	Mikrochirurgische Instrumente / Microsurgical Instruments	114
/	Raspatorien / Raspatories	114
/	Handinstrumente / Manual instruments	115
/	Chirurgieinstrumente / Surgery Instruments	116-127
/	Edelstahl-Bohrerstände / Stainless Steel Bur Blocks	128-129
/	Trays	130

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Meisinger stands for high quality medical devices since 1888. The quality management system of a company which manufactures medical devices must meet specific special requirements. These extremely high requirements are defined in ISO 13485 and meticulously complied with by our company. A MDSAP certificate according to ISO 13485:2016 confirms compliance with the requirements of international authorities in the USA (FDA), Canada (Health Canada), Australia (TGR), Japan (MHLW) and Brazil (ANVISA). All medical devices which you purchase from us as customer, comply with all applicable requirements of the Medical Device Directive 93/42/EEC. Our company is certified by an independent Notified Body and certification is performed according to the specifications of standards. Current certificates can be found on our homepage www.meisinger.de

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All rights reserved.

Reproduction, also by extract and reproducing photos, are only permitted with written authorization of Hager & Meisinger GmbH, Neuss.

All offers, orders, and deliveries are subject to Meisinger's "General Sales and Delivery Terms."

We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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Membranes, Sutures, Pro-Fix™

HERSTELLER MANUFACTURER

Osteogenics Biomedical, Inc.
4620 71st Street | Building 78-79
Lubbock, TX 79424 | USA
Phone: 1.888.796.1923
Fax: 806.796.0059
E-Mail: sales@osteogenics.com
Internet: www.osteogenics.com

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Cytoplast™ and Pro-Fix™ are registered trademarks of Osteogenics Biomedical, Inc..
cerabone® and collafleece® are registered trademarks of botiss biomaterial GmbH. .
NanoBone® ist ein eingetragenes Warenzeichen der Biocomposites GmbH. NanoBone® is a registered trademark of Biocomposites GmbH.
botiss biomaterials GmbH
Biocomposites GmbH
Cytoplast™ und Pro-Fix™ sind eingetragene Warenzeichen der Osteogenics Biomedical, Inc..